

Product Name:	ABT-HC-CP-16
Description:	<p>Premier Glass Door Chromatography Refrigerator  16 Cu. Ft. capacity  One swing glass door, self-closing, right hinged, not reversible  Microprocessor Temperature Controller  Adjustable operating temperature range: 1°-10°C  Digital Temperature Display and Alarm Module  Sample simulated digital temperature display; .1° Resolution  Adjustable temperature display in Celsius or Fahrenheit  Audible and visual high/low temperature alarms  Sensor error alarm  Power failure alarm  Resettable 1-point Min/Max History  Remote alarm contacts  2 probes; 1 in air and 1 in sample bottle  Alarm validation (audible, visual, and remote alarm)  Alarm mute with ring back  Three adjustable shelves (adjustable in ½" increments)  Forced draft air circulation  Keyed door lock  One GFCI protected duplex electrical outlet  Two 2" access ports with spring-loaded, self-closing gasketed metal covers  Magnetic door gasket for positive seal  Heavy duty door handle  LED interior lights are safety shielded and switch controlled  Swivel Casters (3 ½")  Probe access port (3/4")  Advanced hydrocarbon detector  White powder-coated steel interior and exterior  High density urethane foam cabinet insulation  Hydrocarbon, natural refrigerant (R290)  Two year parts and labor warranty, plus an additional three compressor parts warranty  Exterior dimensions: 26 1/2" W* x 29 3/4" D x 79" H with casters; 75 1/2" H without casters. *Width w/access ports closed; 34 1/2" width w/access ports open  Cycle defrost  115V, 60 Hz, 3 Amps, 1/5 HP  ENERGY STAR CERTIFIED  ETL listed  Shipping weight: 283 lbs. Note: This unit requires two (2) outlets</p>