

Product Name:	ABT-HC-CS-36
Description:	<p>Standard Glass Door Chromatography Refrigerator  36 Cu. Ft. capacity  Two swing glass doors, self-closing, right hinged, not reversible  Microprocessor temperature controller  Digital temperature display  Audible and visual high and low temperature alarms  2 probes; 1 in air and 1 in sample bottle  Remote alarm contacts  Adjustable operating temperature range: 1°C – 10°C  Keyed door locks  One GFCI protected duplex electrical outlet  One 2" access port with spring-loaded, self-closing gasketed metal cover  Magnetic door gasket for positive seal  Eight adjustable shelves (adjustable in 1/2" increments)  Forced draft air circulation  Heavy duty door handles  LED interior lights are safety shielded and switch controlled  Swivel Casters (3 1/2")  Probe access port (3/8")  Advanced hydrocarbon detector  White powder-coated steel exterior and interior  High density urethane foam cabinet insulation  Hydrocarbon, natural refrigerant (R290)  One year parts and labor warranty, plus an additional four year compressor parts warranty  Exterior dimensions: 40 3/8" W x 34 3/4" D x 81 3/4" H with casters; 78 1/4" H without casters  Cycle defrost  115V, 60 Hz, 4.5 Amps, 1/4 HP  ETL listed  Shipping weight: 405 lbs. Note: This unit requires two (2) outlets</p>