

Product Name:	ABT-HC-FRFC10
Description:	<p>Standard Flammable Refrigerator &amp; Freezer Combination  9 Cu. Ft. capacity (5 Ref./4 Freezer)  Two exterior solid swing doors right hinged, not reversible  Microprocessor temperature controller, refrigerator and freezer  Digital temperature display, refrigerator and freezer  Adjustable operating temperature range: Refrigerator: 1°C – 10°C; Freezer: -15°C – -25°C  Audible and visual high and low temperature alarms, refrigerator and freezer  2 probes; 1 in air and 1 in sample bottle, in both Refrigerator and Freezer  Remote alarm contacts, refrigerator and freezer  3 adjustable shelves refrigerator/2 fixed shelves (with 3 removable baskets) freezer  Keyed door locks  Magnetic door gasket for positive seal  Recessed door handle, refrigerator and freezer  Leveling legs  Stacking kit  Two probe access ports (3/8"); one, refrigerator/one, freezer  White powder coated exterior  HFC-free refrigerant (R600a), foam insulation and packaging  Two year parts and labor warranty  Exterior dimensions: 23 5/8"W x 24 1/4"D x 65 3/8"H  Manual defrost freezer / cycle defrost refrigerator  115V, 60 Hz, 1.2 Amps Refrigerator and 1.4 Amps Freezer, 1/9 HP Refrigerator and 1/8 HP Freezer  UL/C-UL listed  Units meet NFPA guidelines 45 and 70  Shipping weight: 205 lbs.  *Please note this unit requires 2 outlets.</p>