

Product Name:	ABT-HC-SSP-72G
Description:	<p>Premier Series Stainless Steel Glass Door Refrigerator 72 Cu. Ft. Capacity</p> <p>Three swing glass doors, self-closing, not reversible</p> <p>12 ea. Adjustable Shelves</p> <p>Digital display microprocessor temperature controller</p> <p>Adjustable operating temperature: 1°C to 10°C (34°F to 50°F)</p> <p>Audible and visual Hi/Lo temperature alarm</p> <p>2 probes; 1 in air and 1 in sample bottle</p> <p>Sensor error alarm</p> <p>Continuous product (simulated) temperature display</p> <p>Remote alarm contacts</p> <p>Forced draft air circulation</p> <p>Keyed door locks</p> <p>Magnetic door gaskets for positive seal</p> <p>Heavy duty door handles</p> <p>LED interior lights are safety shielded and switch controlled</p> <p>Probe access port (3/4")</p> <p>Stainless steel exterior</p> <p>White-powder coated aluminum interior with 304 stainless steel floor</p> <p>Easily cleaned smooth interior</p> <p>Hydrocarbon natural refrigerant (R290)</p> <p>2 year parts and labor warranty, plus an additional 3 year compressor parts warranty</p> <p>6 swivel casters (5 3/4")</p> <p>Exterior dimensions: 81" W x 35 3/4" D x 83 1/4" H with casters; 77 1/2" H without casters</p> <p>Cycle defrost</p> <p>115V, 60 Hz, 7.5 Amps, 1/2 HP</p> <p>ETL listed</p> <p>Shipping weight: 786 lbs.</p>