

Product Name:	ABT-HC-UCFS-0504W
Description:	<p>Standard Undercounter Refrigerator Freestanding 5 Cu. Ft. Capacity One solid swing door, right hinged, not reversible Mechanical thermostat Temperature range : 0°C to 10°C Three shelves (1 fixed/2 adjustable) Keyed door lock Magnetic door gasket for positive seal Recessed door handle Leveling legs Interior light Probe access port (3/8") White cabinet and door High density urethane foam cabinet and door insulation HFC-free refrigerant (R600a), foam insulation and packaging One year parts and labor warranty Exterior dimensions: 20 1/4"W x 22" D x 31 1/2"H <i>Note: This unit must have 4" clearance on sides and back for adequate ventilation</i> Cycle defrost 115V, 60 Hz, .90 Amps, 1/7 HP UL/C-UL listed Shipping weight: 86 lbs.</p>