

Product Name:	AMIV-1 PS
Description:	<p>ABS Auto Max System with Polycarbonate Package System Auto Fill Tank Total Vial Capacity (2 mL vials): 10,400 Total Bag Capacity: 235 LN2 Capacity (Liters): 165 Static Evaporation Rate (Liters/Day): 5 Normal working time (days): 33 Power requirement: 24 Vac Exterior Width (mm): 586.7 Overall Height (mm): 1117.6 Exterior Depth (mm): 774.7 Usable Height (mm): 736.6 Internal Diameter (mm): 533.4 Sentry Controller</p> <p>Controller features:</p> <ul style="list-style-type: none"> ○ Temperature display (Celsius) ○ Level accuracy with 4 direct read sensors ○ Large format LED display ○ Alarm ring back ○ Redundant level sensor system ○ Remote monitoring system available ○ Liquid level indicator ○ Alarm package (audible, visual, and remote) including level (high and low), high temperature, and sensor failure <p>UL-Listed transformer Stainless steel vacuum tank Recessed casters Fiberglass top Full access opening Leveling platform Stainless steel armor-plated transfer hose 7-inch foam insulated hinged lid with locking provision, protective boot, and pneumatic lifters One-year warranty on all parts and labor Three-year vacuum warranty</p> <p>Package System Includes:</p> <ul style="list-style-type: none"> - (7) 13-2 Tower Racks (81 or 100 cell boxes) - (4) 13-2 Mini Tower Racks (25 cell mini boxes) - (91) 81 or 100 Cell Polycarbonate Boxes - (54) 25 Cell Mini Polycarbonate Boxes - (1) LN₂ Measuring Ruler - (1) Pair CryoGloves