

Product Name:	PH-ABT-HC-RFC1020G
Description:	Premier Pharmacy Refrigerator & Freezer Combination (-20°C) Freestanding 5.2 cu. ft. capacity (Refrigerator) / 4 cu. ft. capacity (Freezer) One glass swing door (fridge), one solid swing door (freezer) Microprocessor temperature controller Digital temperature display Audible and visual high and low temperature alarms 2 probes; 1 in air and 1 in sample bottle Remote alarm contacts Adjustable operating temperature range - Freezer: 5°F – -13°F (-15° C – -25°C); Refrigerator: 36°F to 46°F (2°C to 8°C) Temperature monitor device (TMD) complies with the current CDC guidelines, with 3 years certification of calibration, “buffered” probe in the product simulated solution, min/max memory. F/C switchable, field installable, and visual & audible temp alarm 2 fixed shelves (freezer); 3 adjustable shelves (fridge) Bottom footed shelf to allow air circulation Back wall evaporator cover Keyed door locks Magnetic door gasket for positive seal Recessed door handles Stacking kit Leveling legs LED interior light Probe access port (3/8") White cabinet and door High density urethane foam cabinet and door insulation HFC-free refrigerant (R600a), foam insulation and packaging Two year parts and labor warranty, plus an additional three year compressor parts only warranty Exterior dimensions: 23 3/4" W x 24"D x 65"H <i>Note: This unit must have 4" clearance on sides and back for adequate ventilation</i> Manual defrost (freezer); Cycle defrost (refrigerator) 115V, 60 Hz 1.4 Amps, 1/8 HP (freezer); 1.3 Amps, 1/6 HP (refrigerator) UL/C-UL listed Shipping weight: 266 lbs.