

Product Name:	PH-ABT-HC-RFC1030G
Description:	<p>Premier Pharmacy Refrigerator &amp; Freezer Combination (-30°C)  Freestanding  5.2 cu. ft. capacity (Refrigerator) / 4 cu. ft. capacity (Freezer) One glass swing door (fridge), one solid swing door (freezer)  Microprocessor temperature controller  Digital temperature display  Audible and visual high and low temperature alarms  2 probes; 1 in air and 1 in sample bottle  Remote alarm contacts  Adjustable operating temperature range:  Freezer: -4°F to -22°F (-20°C to -30°C). Factory Preset at -30°C;  Refrigerator: 36°F to 46°F (2°C to 8°C)  Temperature monitor device (TMD) complies with the current CDC guidelines, with 3 years certification of calibration, "buffered" probe in the product simulated solution, min/max memory. F/C switchable, field installable, and visual &amp; audible temp alarm  3 adjustable shelves (freezer); 3 adjustable shelves (fridge)  Bottom footed shelf to allow air circulation  Back wall evaporator cover  Keyed door locks  Magnetic door gasket for positive seal  Recessed door handles  Stacking kit  Leveling legs  LED interior light  Probe access port (3/8")  White cabinet and door  High density urethane foam cabinet and door insulation  HFC-free refrigerant (R600a) Refrigerator; Natural, Hydrocarbon refrigerant (R290) Freezer  Two year parts and labor warranty, plus an additional three year compressor parts only warranty  Exterior dimensions: 23 3/4" W x 24"D x 65"H  <i>Note: This unit must have 4" clearance on sides and back for adequate ventilation</i>  Manual defrost (freezer); Cycle defrost (refrigerator)  115V, 60 Hz  2.2 Amps, 1/4 HP (freezer); 1.3 Amps, 1/6 HP (refrigerator)  UL/C-UL listed  Shipping weight: 268 lbs.</p>