

Product Name:	PH-ABT-HC-S16S
Description:	<p>Standard Pharmacy/Vaccine Solid Door Refrigerator  16 Cu. Ft. capacity  One swing solid door, self-closing, right hinged, not reversible  Microprocessor temperature controller  Digital temperature display  Audible and visual high and low temperature alarms  Remote alarm contacts  2 probes; 1 in air and 1 in sample bottle  Adjustable operating temperature range: 36°F – 46°F (2°C – 8°C)  Temperature monitor device (TMD) complies with the current CDC guidelines, with 3 years certification of calibration, “buffered” probe in the product simulated solution, min/max memory. F/C switchable, field installable, and visual &amp; audible temp alarm  Three adjustable shelves (adjustable in ½” increments)  Vaccine storage power cord warning label  Forced draft air circulation  Heavy duty door handle  Keyed door lock  Magnetic door gasket for positive seal  LED interior lights are safety shielded and switch controlled  Swivel Casters (3 ½”)  Probe access port (⅜”)  White powder-coated steel exterior and interior  Hydrocarbon, natural refrigerant (R290)  Pharmacy refrigerator/freezer toolkit and temperature logs  One year parts and labor warranty, plus an additional four year compressor parts warranty  Exterior dimensions: 25” W x 29 ¾” D x 79” H with casters; 75 ½” H without casters  Cycle defrost  115V, 60 Hz, 3 Amps, 1/5 HP  ENERGY STAR CERTIFIED  ETL listed  Shipping weight: 277 lbs.</p>