

Product Name:	PH-ABT-HC-UCBI-0420SS
Description:	Premier Pharmacy/Vaccine Undercounter Freezer Built-In 4.2 Cu. Ft. Capacity One stainless steel swing door, right hinged* Microprocessor temperature controller Digital temperature display Audible and visual high and low temperature alarms Remote alarm contacts Adjustable operating temperature range: 5°F – -13°F (-15°C – -25°C) Temperature monitor device (TMD) complies with the current CDC guidelines, with 3 years certification of calibration, “buffered” probe in the product simulated solution, min/max memory. F/C switchable, field installable, and visual & audible temp alarm Two fixed shelves Vaccine storage power cord warning label Keyed door lock Recessed door handle Magnetic door gasket for positive seal Leveling legs Probe access port (3/8") Stainless steel exterior High density urethane foam cabinet and door insulation HFC-free refrigerant (R600a), foam insulation and packaging Pharmacy refrigerator/freezer toolkit and temperature logs Two year parts and labor warranty, plus an additional three year compressor parts only warranty Pyxis®, Omnicell® and AcuDose RX® compatible Exterior dimensions: 23 ¾ "W x 24 ½" D x 33 3/8 "H Manual defrost 115V, 60 Hz, 1.5 Amps, 1/8 HP Energy Consumption 0.97 kWh/day UL/C-UL Listed ENERGY STAR CERTIFIED Shipping weight: 131 lbs. *This product is available left hinged, product number PH-ABT-HC-UCBI-0420SS-LH