

Product Name:	PH-ABT-HC-UCFS-0120-LH
Description:	<p>Premier Pharmacy/Vaccine Countertop Laboratory Freezer  Freestanding-LEFT HINGED  1.7 Cu. Ft. Capacity  One solid swing door, left hinged, not reversible  Microprocessor temperature controller  Digital temperature display  Audible and visual high and low temperature alarms  Remote alarm contacts  2 probes; 1 in air and 1 in glycol bottle  Adjustable operating temperature range: 5°F – -13°F (-15°C – -25°C)  Temperature monitor device (TMD) complies with the current CDC guidelines, with 3 years certification of calibration, “buffered” probe in the product simulated solution, min/max memory, F/C switchable, field installable, and visual &amp; audible temp alarm  2 adjustable shelves with guard rail to prevent product from being pushed against back wall  Bottom footed shelf to allow air circulation  Vaccine storage power cord warning label  Keyed door lock  Magnetic door gasket for positive seal  Leveling legs  Integrated door handle  Probe access port (3/8")  Powder coat exterior  High density urethane foam cabinet and door insulation  HFC-free refrigerant (R600a), foam insulation and packaging  Pharmacy refrigerator/freezer toolkit and temperature logs  Two year parts and labor warranty, plus an additional three year compressor parts only warranty  Pyxis®, Omnicell® and AcuDose RX® compatible  Exterior dimensions: 23 ¾"W x 24"D x 20 ⅜"H  <i>Note: This unit must have 4" clearance on sides and back for adequate ventilation</i>  Manual defrost  115V, 60 Hz, 1 Amp, 1/8 HP  Energy Consumption 0.75 kWh/day  UL/C-UL listed  ENERGY STAR CERTIFIED  Shipping weight: 106 lbs.</p>