| | |
|---|---|
| Product Name: | ABT-HC-10PG |
| | |
| Description: | Premier Glass Door Compact Laboratory Refrigerator<br>10.5 Cu. Ft. Capacity<br>One glass swing door, right hinged, not reversible<br>Microprocessor temperature controller<br>Digital temperature display with sample simulated probe<br>Audible and visual high and low temperature alarms<br>Remote alarm contacts<br>Adjustable operating temperature range: 1°C to 10°C<br>Six adjustable shelves with guard rail on back to prevent product from being pushed against back wall<br>Bottom footed shelf to allow air circulation<br>Forced draft air circulation<br>Back wall evaporator cover<br>Keyed door lock<br>Magnetic door gasket for positive seal<br>Leveling legs<br>LED interior light<br>Probe access port (3/8")<br>Powder coat exterior<br>High density urethane foam cabinet insulation<br>HFC-free refrigerant (R600a), foam insulation and packaging<br>Two year parts and labor warranty, plus an additional three year compressor parts warranty<br>Exterior dimensions: 23 ¾"W x 26 ½"D x 59"H<br>*Note: This unit must have 4" clearance on sides and back for adequate ventilation*<br>UL/C-UL listed<br>Cycle defrost<br>115V, 60 Hz, 1/3 HP, 2.3 Amps<br>Energy Consumption 2.42 kWh/day<br>Shipping weight: 209 lbs. |